

3 Tesla Users Primer and Training Sessions

Karen E. Bove Bettis, RT (R)(MR)

MRI Technologist and Research Assistant

Paula Rowser, RT (R)(MR)

MRI Technologist and Research Assistant

Objectives

- Review the guidelines for human scanning and supply references for further information
- To scan safely
- To scan efficiently
- To find the scan data

Lesson 3

Data and Image Transfer

Lesson 3

Administrative Issues and Procedures

Safety Issues

Scan Operations and Procedures


- Data and Image Transfer
 - EPI Data
 - Raw Files (P Files)
 - GE Real Time Files
 - Image Files
 - GINX
 - Signa Archive

Scan Room Peripheral Equipment

Scan Room Issues

Miscellaneous

Advanced Training






Data and Image Transfers

- Scanner Console
 - Image Files
 - (Browser)
 - P File Numbers
 - (Tools_C Shell)
- GE Sun ADW
 - GE Real Time
 - GINX
- Potassium
 - GE Real Time
 - P Files



The GE Sun ADW on 3T_1

“Lay of the Land”




P Files or GE Real Time Data

- GINX
 - xferx
 - Image Data
- NIH EPI
 - P Files
 - sdc-adw
 - Alias = raw
- Image Data
 - Anatomical data
 - Structural data
 - Clinical data
- GE RealTime
 - sdc-nfs
 - Alias = epirt

Signa LX Console Image Files

Scanner Operations Icon



- Date and Time
- Informix space
- Message Board
- **Scanner Operations**
 - Patient Information
 - Protocol RX
- Protocol Manager
- Browser
- Network
- Tools
- Piece of the Puzzle

LX Console > Image Files

Image scanning =
non EPI or SPEP


- Anatomical scans
- Structural scans
- Clinical Scans



What Are Image Files?

For example:

- 2d & 3d FSPGR
- 2d & 3d SPGR
- SE
- FSE
- FLAIR
- MRA
- MP-RAGE
- etc



The Browser Icon



- Date and Time
 - Hard Drive Space
 - Message Board
 - Scanner Operations
 - Protocol Manager
-
- Browser
 - Archive functions
 - Filming functions
 - Network (local) functions
-
- Network
 - Tools
 - Piece of the Puzzle

Save Image Files to Archive Media

Application Selection Remove Sort Network Archive Queue Utilities Services Messages

Examinations :

Exam	Name	Date	Description	Mod	Fmt	A
7254/	3T crotty, laura	Jul 14 02	tools_tellin	MR	Adv	Y
7253/	3T jenkins, tereza	Jul 14 02	tools_tellin	MR	Adv	Y
7252/	3T jenkins, tereza	Jul 14 02	tools_tellin	MR	Adv	Y
7251/	3T crawford, andrew	Jul 14 02	tools_tellin	MR	Adv	Y
7250/	3T crawford, andrew	Jul 14 02	tools_tellin	MR	Adv	Y
02	FMRI brain	MR	Adv	Y		
02	FMRI brain	MR	Adv	Y		
02	FMRI brain	MR	Adv	Y		
02	FMRI brain	MR	Adv	Y		
02	Morphometry	MR	Adv	Y		
02	FMRI lang	MR	Adv	Y		
02	fmri languag	MR	Adv	Y		
02	FMRI lang	MR	Adv	Y		
02	Moprhometry	MR	Adv	Y		
92	SMPTE	MR	Adv	N		

Exam no 1000, Jul 21 92, SMPTE

Ser	Type	Imgs	Description	Mod	Manf
1	SSAVE	1		MR	GEMS

January 14 7:27 23730 left @ 256² 5921 left @ 512²

Disk 24% full

Image quality may be degraded. Please see message log.

Browser

Idle

5666/1/3 3/3

Media attached

Printed Ex: 5652 Se: 4 Im: 91

Network Status

GE

Image Files

- ALL image files should be saved to the Signa MOD
 - Click on Browser
 - Click on Patient
 - Click Archive
 - Click Save Exam




Image files can be restored to Signa and re-Ginxed at later dates **IF** saved to MOD first. All image data deleted before saved MOD to the is lost forever!

Image Files

- Transfer the image files needed to the Sun ADW via the GINX procedure
- Use rmb
 - IDBM Apps
 - GINX
 - click on the patient
 - Click *Get from Remote*



Signa LX Console P Files

What are P Files ?

**** Raw data plus the header information ****

- Save raw data for recon at a later date
 - Researcher specific
 - Venograms
 - Technical testing
- Protocols that use P Files for Reconstruction
 - NIH EPI
 - SPEP
 - EPI j bEPI

The Tools Icon



- Date and Time
- Informix space
- Message Board
- Scanner Operations
- Protocol Manager
- Browser
- Network
- Tools
 - Engineer Options
 - Reset TPS
 - C Shell
 - System Shutdown
- Piece of the Puzzle

Scanner Console_P Files



- Click on the Tools Icon

Scanner Console_P Files



- Click on the Tools Icon
- Click on Cshell

Scanner Console_P Files



- Click on the Tools Icon
- Click on Cshell

At Prompt:


/usr/g/bin is the default)

Type:

`cd /usr/g/mrraw`

- Hit enter

P Files – Change Directory


A terminal window titled 'winterm' with a dark blue background and white text. It shows two shell prompts. The first prompt is '{sdc@fim3T}[1] ' followed by the command 'pwd'. The output is '/usr/g/bin'. The second prompt is '{sdc@fim3T}[2] ' followed by the command 'cd /usr/g/mrraw'. A green cursor is visible at the end of the second command.

```
winterm
{sdc@fim3T}[1] pwd
/usr/g/bin
{sdc@fim3T}[2] cd /usr/g/mrraw
```

cd into:

`/usr/g/mrraw`

Scanner Console_P Files



```
winterm
{sdC@fim3T}[1] pwd
/usr/g/bin
{sdC@fim3T}[2] cd /usr/g/mrraw
```

`cd /usr/g/mrraw`


Hit enter

Type:

`ls -lrt P*`

- Hit enter

P Files Typing > ls -lrt **VS** ls -lrt P*



```
winterm
{sdc@fim3T}[1] pwd
/usr/g/bin
{sdc@fim3T}[2] cd /usr/g/mrraw
{sdc@fim3T}[3] ls -lrt
```

- ls -lrt

P Files List_Typing ls -lrt

```
winterm
lrwxrwxr-x 1 sdc informix 8 Jul 5 06:38 01030 -> P01024.7*
lrwxrwxr-x 1 sdc informix 8 Jul 5 06:38 01031 -> P01024.7*
-rwxrwxr-x 1 sdc informix 71754800 Jul 5 06:51 P02560.7*
lrwxrwxr-x 1 sdc informix 8 Jul 5 06:51 02560 -> P02560.7*
lrwxrwxr-x 1 sdc informix 8 Jul 5 06:51 02561 -> P02560.7*
lrwxrwxr-x 1 sdc informix 8 Jul 5 06:51 02562 -> P02560.7*
lrwxrwxr-x 1 sdc informix 8 Jul 5 06:51 02563 -> P02560.7*
lrwxrwxr-x 1 sdc informix 8 Jul 5 06:51 02564 -> P02560.7*
lrwxrwxr-x 1 sdc informix 8 Jul 5 06:51 02565 -> P02560.7*
lrwxrwxr-x 1 sdc informix 8 Jul 5 06:51 02566 -> P02560.7*
lrwxrwxr-x 1 sdc informix 8 Jul 5 06:51 02567 -> P02560.7*
-rwxrwxr-x 1 sdc informix 1398832 Jul 5 07:14 P04096.7*
lrwxrwxr-x 1 sdc informix 8 Jul 5 07:14 04096 -> P04096.7*
-rwxrwxr-x 1 sdc informix 1398832 Jul 5 07:18 P05120.7*
lrwxrwxr-x 1 sdc informix 8 Jul 5 07:18 05120 -> P05120.7*
-rwxrwxr-x 1 sdc informix 1398832 Jul 5 07:23 P06144.7*
lrwxrwxr-x 1 sdc informix 8 Jul 5 07:23 06144 -> P06144.7*
-rwxrwxr-x 1 sdc informix 13625904 Jul 5 07:32 P07168.7*
lrwxrwxr-x 1 sdc informix 8 Jul 5 07:32 07168 -> P07168.7*
lrwxrwxr-x 1 sdc informix 8 Jul 5 07:32 07169 -> P07168.7*
lrwxrwxr-x 1 sdc informix 8 Jul 5 07:32 07170 -> P07168.7*
lrwxrwxr-x 1 sdc informix 8 Jul 5 07:32 07171 -> P07168.7*
lrwxrwxr-x 1 sdc informix 8 Jul 5 07:32 07172 -> P07168.7*
lrwxrwxr-x 1 sdc informix 8 Jul 5 07:32 07173 -> P07168.7*
lrwxrwxr-x 1 sdc informix 8 Jul 5 07:32 07174 -> P07168.7*
lrwxrwxr-x 1 sdc informix 8 Jul 5 07:32 07175 -> P07168.7*
lrwxrwxr-x 1 sdc informix 8 Jul 5 07:32 07176 -> P07168.7*
lrwxrwxr-x 1 sdc informix 8 Jul 5 07:32 07177 -> P07168.7*
-rwxrwxr-x 1 sdc informix 1398832 Jul 5 07:36 P08704.7*
lrwxrwxr-x 1 sdc informix 8 Jul 5 07:36 08704 -> P08704.7*
-rwxrwxrwx 1 sdc informix 320000 Jul 5 08:10 HEADER_POOL*
{sdc@fim3T}[4]
```

A confusing order results Results when typing:

ls lrt

P Files – Date and Time Stamp

```
winterm
-rwxrwxr-x 1 sdc informix 1744432 Jul 3 19:29 P38912.7*
-rwxrwxr-x 1 sdc informix 24366640 Jul 3 19:38 P39424.7*
-rwxrwxr-x 1 sdc informix 24366640 Jul 3 19:44 P39936.7*
-rwxrwxr-x 1 sdc informix 24366640 Jul 3 19:56 P40960.7*
-rwxrwxr-x 1 sdc informix 24366640 Jul 3 20:01 P41472.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 21:49 P50688.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 21:51 P51200.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 21:52 P51712.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 23:41 P59904.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 23:42 P60416.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 23:43 P60928.7*
-rwxrwxr-x 1 sdc informix 1877040 Jul 4 13:48 P04608.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 14:00 P05632.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 14:23 P06656.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 14:41 P07680.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 14:52 P08192.7*
-rwxrwxr-x 1 sdc informix 1877040 Jul 4 15:18 P11776.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 15:29 P12288.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 15:40 P12800.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 15:52 P13312.7*
-rwxrwxr-x 1 sdc informix 42574896 Jul 4 15:57 P13824.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 16:09 P14336.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 16:21 P14848.7*
-rwxrwxr-x 1 sdc informix 42574896 Jul 4 16:31 P16896.7*
-rwxrwxr-x 1 sdc informix 71754800 Jul 5 06:38 P01024.7*
-rwxrwxr-x 1 sdc informix 71754800 Jul 5 06:51 P02560.7*
-rwxrwxr-x 1 sdc informix 1398832 Jul 5 07:14 P04096.7*
-rwxrwxr-x 1 sdc informix 1398832 Jul 5 07:18 P05120.7*
-rwxrwxr-x 1 sdc informix 1398832 Jul 5 07:23 P06144.7*
-rwxrwxr-x 1 sdc informix 13625904 Jul 5 07:32 P07168.7*
-rwxrwxr-x 1 sdc informix 1398832 Jul 5 07:36 P08704.7*
{sdc@fim3T}[4]
```

Orderly results when typing:

```
ls -lrt P*
```


C Shell_(ls -lrt P*)

```
winterm
-rwxrwxr-x 1 sdc informix 1744432 Jul 3 19:29 P38912.7*
-rwxrwxr-x 1 sdc informix 24366640 Jul 3 19:38 P39424.7*
-rwxrwxr-x 1 sdc informix 24366640 Jul 3 19:44 P39936.7*
-rwxrwxr-x 1 sdc informix 24366640 Jul 3 19:56 P40960.7*
-rwxrwxr-x 1 sdc informix 24366640 Jul 3 20:01 P41472.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 21:49 P50688.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 21:51 P51200.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 21:52 P51712.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 23:41 P59904.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 23:42 P60416.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 23:43 P60928.7*
-rwxrwxr-x 1 sdc informix 1877040 Jul 4 13:48 P04608.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 14:00 P05632.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 14:23 P06656.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 14:41 P07680.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 14:52 P08192.7*
-rwxrwxr-x 1 sdc informix 1877040 Jul 4 15:18 P11776.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 15:29 P12288.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 15:40 P12800.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 15:52 P13312.7*
-rwxrwxr-x 1 sdc informix 42574896 Jul 4 15:57 P13824.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 16:09 P14336.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 16:21 P14848.7*
-rwxrwxr-x 1 sdc informix 42574896 Jul 4 16:31 P16896.7*
-rwxrwxr-x 1 sdc informix 71754800 Jul 5 06:38 P01024.7*
-rwxrwxr-x 1 sdc informix 71754800 Jul 5 06:51 P02560.7*
-rwxrwxr-x 1 sdc informix 1398832 Jul 5 07:14 P04096.7*
-rwxrwxr-x 1 sdc informix 1398832 Jul 5 07:18 P05120.7*
-rwxrwxr-x 1 sdc informix 1398832 Jul 5 07:23 P06144.7*
-rwxrwxr-x 1 sdc informix 13625904 Jul 5 07:32 P07168.7*
-rwxrwxr-x 1 sdc informix 1398832 Jul 5 07:36 P08704.7*
{sdc@fim3T}[4]
```

File Size

C Shell_(ls -lrt P*)

```
winterm
-rwxrwxr-x 1 sdc informix 1744432 Jul 3 19:29 P38912.7*
-rwxrwxr-x 1 sdc informix 24366640 Jul 3 19:38 P39424.7*
-rwxrwxr-x 1 sdc informix 24366640 Jul 3 19:44 P39936.7*
-rwxrwxr-x 1 sdc informix 24366640 Jul 3 19:56 P40960.7*
-rwxrwxr-x 1 sdc informix 24366640 Jul 3 20:01 P41472.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 21:49 P50688.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 21:51 P51200.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 21:52 P51712.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 23:41 P59904.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 23:42 P60416.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 23:43 P60928.7*
-rwxrwxr-x 1 sdc informix 1877040 Jul 4 13:48 P04608.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 14:00 P05632.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 14:23 P06656.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 14:41 P07680.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 14:52 P08192.7*
-rwxrwxr-x 1 sdc informix 1877040 Jul 4 15:18 P11776.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 15:29 P12288.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 15:40 P12800.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 15:52 P13312.7*
-rwxrwxr-x 1 sdc informix 42574896 Jul 4 15:57 P13824.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 16:09 P14336.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 16:21 P14848.7*
-rwxrwxr-x 1 sdc informix 42574896 Jul 4 16:31 P16896.7*
-rwxrwxr-x 1 sdc informix 71754800 Jul 5 06:38 P01024.7*
-rwxrwxr-x 1 sdc informix 71754800 Jul 5 06:51 P02560.7*
-rwxrwxr-x 1 sdc informix 1398832 Jul 5 07:14 P04096.7*
-rwxrwxr-x 1 sdc informix 1398832 Jul 5 07:18 P05120.7*
-rwxrwxr-x 1 sdc informix 1398832 Jul 5 07:23 P06144.7*
-rwxrwxr-x 1 sdc informix 13625904 Jul 5 07:32 P07168.7*
-rwxrwxr-x 1 sdc informix 1398832 Jul 5 07:36 P08704.7*
{ sdc@fim3T } [4]
```

Date Acquired

C Shell_(ls -lrt P*)

```
winterm
-rwxrwxr-x 1 sdc informix 1744432 Jul 3 19:29 P38912.7*
-rwxrwxr-x 1 sdc informix 24366640 Jul 3 19:38 P39424.7*
-rwxrwxr-x 1 sdc informix 24366640 Jul 3 19:44 P39936.7*
-rwxrwxr-x 1 sdc informix 24366640 Jul 3 19:56 P40960.7*
-rwxrwxr-x 1 sdc informix 24366640 Jul 3 20:01 P41472.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 21:49 P50688.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 21:51 P51200.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 21:52 P51712.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 23:41 P59904.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 23:42 P60416.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 23:43 P60928.7*
-rwxrwxr-x 1 sdc informix 1877040 Jul 4 13:48 P04608.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 14:00 P05632.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 14:23 P06656.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 14:41 P07680.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 14:52 P08192.7*
-rwxrwxr-x 1 sdc informix 1877040 Jul 4 15:18 P11776.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 15:29 P12288.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 15:40 P12800.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 15:52 P13312.7*
-rwxrwxr-x 1 sdc informix 42574896 Jul 4 15:57 P13824.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 16:09 P14336.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 16:21 P14848.7*
-rwxrwxr-x 1 sdc informix 42574896 Jul 4 16:31 P16896.7*
-rwxrwxr-x 1 sdc informix 71754800 Jul 5 06:38 P01024.7*
-rwxrwxr-x 1 sdc informix 71754800 Jul 5 06:51 P02560.7*
-rwxrwxr-x 1 sdc informix 1398832 Jul 5 07:14 P04096.7*
-rwxrwxr-x 1 sdc informix 1398832 Jul 5 07:18 P05120.7*
-rwxrwxr-x 1 sdc informix 1398832 Jul 5 07:23 P06144.7*
-rwxrwxr-x 1 sdc informix 13625904 Jul 5 07:32 P07168.7*
-rwxrwxr-x 1 sdc informix 1398832 Jul 5 07:36 P08704.7*
{ sdc@fim3T } [4]
```

Time
Acquired


C Shell_(ls -lrt P*)

```
winterm
-rwxrwxr-x 1 sdc informix 1744432 Jul 3 19:29 P38912.7*
-rwxrwxr-x 1 sdc informix 24366640 Jul 3 19:38 P39424.7*
-rwxrwxr-x 1 sdc informix 24366640 Jul 3 19:44 P39936.7*
-rwxrwxr-x 1 sdc informix 24366640 Jul 3 19:56 P40960.7*
-rwxrwxr-x 1 sdc informix 24366640 Jul 3 20:01 P41472.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 21:49 P50688.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 21:51 P51200.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 21:52 P51712.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 23:41 P59904.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 23:42 P60416.7*
-rwxrwxr-x 1 sdc informix 2630704 Jul 3 23:43 P60928.7*
-rwxrwxr-x 1 sdc informix 1877040 Jul 4 13:48 P04608.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 14:00 P05632.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 14:23 P06656.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 14:41 P07680.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 14:52 P08192.7*
-rwxrwxr-x 1 sdc informix 1877040 Jul 4 15:18 P11776.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 15:29 P12288.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 15:40 P12800.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 15:52 P13312.7*
-rwxrwxr-x 1 sdc informix 42574896 Jul 4 15:57 P13824.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 16:09 P14336.7*
-rwxrwxr-x 1 sdc informix 84968496 Jul 4 16:21 P14848.7*
-rwxrwxr-x 1 sdc informix 42574896 Jul 4 16:31 P16896.7*
-rwxrwxr-x 1 sdc informix 71754800 Jul 5 06:38 P01024.7*
-rwxrwxr-x 1 sdc informix 71754800 Jul 5 06:51 P02560.7*
-rwxrwxr-x 1 sdc informix 1398832 Jul 5 07:14 P04096.7*
-rwxrwxr-x 1 sdc informix 1398832 Jul 5 07:18 P05120.7*
-rwxrwxr-x 1 sdc informix 1398832 Jul 5 07:23 P06144.7*
-rwxrwxr-x 1 sdc informix 13625904 Jul 5 07:32 P07168.7*
-rwxrwxr-x 1 sdc informix 1398832 Jul 5 07:36 P08704.7*
{sdc@fim3T}[4]
```

P File
Numerical
Name

Using NIH EPI

- Log on to Potassium
 - sdc-adw
 - pidof demon2
- Open terminal on Sun ADW
 - xhost +
- Scan on Signa using:
 - jbEPI
 - SPEP
 - NIH 3



AFNI Pop-Up Window
on the ADW

Scanning with NIH EPI

Recording the File Names and Seeing the Images

On Signa:

- Click on the tools icon
 - Open a Cshell
 - cd to:
`/usr/g/mrraw`
 - `ls -lrt P*`
- It is important to write down and to keep track of the P File Names (Numerical) after each run.
- Remember to logout of Potassium by clicking on **fmri done**.

NIH EPI_Potassium

With the following sequences:

NIH EPI

jbEPi

SPEP

The script is automatically making the brik files and sending them to AFNI

NIH EPI_Potassium

When finished scanning, go to Potassium and click on **fmri done** button.

This will automatically make the tar file

Type `ls-lrt` in the window

Note the file name (Jul_25_1340)

Logout of Potassium

Finding the P Files on Potassium

- Reside on sdc-adw

- Type raw (alias)

- ls -lrt

- cd Jun_30_XXX

- P file #
 - FILE_LIST

The P Files will remain on Potassium for 5-7 days only.

It is the researcher's responsibility to copy the files to their own computer.

For example: Jun_30_085

Sun ADW

The Sun ADW

- Connected to the Signa IDBM Database
- Resides to the right of the scan console
- AFNI Pop-up with NIH Epi
- GINX



SUN ADW Screen

fim3t-adw [102]%

Xhost +

Access control disabled, clients can connect from any host



console

/bin/csh

/bin/csh

NIH EPI _ Sun ADW

- **On the GE Sun ADW**
- Open a Terminal window
- Type:
 xhost +
- In this manner, AFNI will automatically pop up on the ADW. Image runs will be listed by the p file name (01024.7).
- When finished scanning, go to Potassium and click on **fmri done** button.
 - Type ls-lrt in the window
 - Note the file name (Jul_25_1340)
 - Logout of Potassium
 - *Transfer the file to personal computer*

Running GE Real Time EPI

- After download of protocol and checking the off-set for fat sat. On the Sun ADW :
- click on:
 - RealTime Control (rmb)
 - Click on Prep scan
- When it fails, click on Manual Prescan (Signa console) and check the prescan parameters
- Write down the values
- Click on Reference Scan (Signa Console)
- Click on Scan (Sun ADW)

SUN ADW Screen

fim3t-adw [102]%

Using the *rmb* click on:

Applications

RealTime Control
Neuro Apps
IDBM Apps
Terminals
Tools
Window Manager
Service Tools



console

/bin/csh

/bin/csh

Tar and Copy Data From ADW

Immediately after scan and before the next patient:

- In a terminal window:
 - Copy (cp) data from SUN through the Network (nfs) to Potassium
- Type: `cd $RTime` and enter
 - `ls -lrt` and enter
 - `cp -rp 36_70_46_1/$arena2/tmp`

(tmp resides on Potassium in sdc-adw)

SUN ADW Screen

fim3t-adw [102]%

Using the *rmb* click on:

Applications

RealTime Control
Neuro Apps
IDBM Apps
Terminals
Tools
Window Manager
Service Tools

IDBM Apps

GINX
Db Maintenance
Clean Failed Image



console

/bin/csh

/bin/csh

SUN ADW

Screen

fim3t-adw [102]%

Using the *rmb* click on:



console

/bin/csh

/bin/csh

Applications

RealTime Control
Neuro Apps
IDBM Apps
Terminals →
Tools
Window Manager
Service Tools

Terminals

Command Tool
Xterm
ShellTool

Colors

White
Purple
Yellow

SUN ADW Screen

fim3t-adw [102]%

Using the *rmb* Click on:

Applications

RealTime Control
Neuro Apps
IDBM Apps
Terminals
Tools
Window Manager
Service Tools

Tools

Text Editor
Calculator
File Manager
Tape Tool
Perfmer
Manual Page Viewer



console

/bin/csh

/bin/csh

SUN ADW Screen

fim3t-adw [102]%

Using the *rmb* Click on:

Applications

RealTime Control
Neuro Apps
IDBM Apps
Terminals
Tools
Window Manager
Service Tools

Window Manager

Redisplay All
Properties
Reload Resource Database
Restart
Lock Screen



console


/bin/csh

/bin/csh


SUN ADW Screen

fim3t-adw [102]%


Using the *rmb* click on:



console



/bin/csh



/bin/csh

Applications

RealTime Control
Neuro Apps
IDBM Apps
Terminals
Tools
Window Manager
Service Tools

Service Tools

Install License Key
Logout

Potassium

Potassium Screen (sdc-adw)

fMRI done

```
Sdc-adw@potassium.nimh.nih.gov:/home/sdc-adw-konsole
```

```
$ pidof demon2
```

```
Sdc-adw@potassium.nimh.nih.gov:/home/sdc-adw-konsole
```

NIH EPI_Potassium

- **On Potassium:**
- At RedHat Login:
 - Login: sdc-adw
 - Password: xxxxxx
- In window:
 - Type:
 - pidof demon2 (should be a number)
- When finished scanning, go to Potassium and click on **fmri done** button.
 - Type ls-lrt in the window
 - Note the file name (Jul_25_1340)
 - Logout of Potassium

Tar Data From Backup (Potassium)

- Telnet to Potassium
 - Login: sdc-nfs
 - Password: xxxxxx
 - Type: epirt and enter
 - Type: ls -lrt and enter
 - Cd into the day and data
 - e.g. > cd DATA_Jul_20/MROR/07301
 - Type: tar cvf 07301.tar 07301/ and enter
 - Type: ls -lrt
 - ftp to own computer

Finding Data After a TPS Reset



- When the TPS is reset on the scanner...
- All the P Files from the given date are transferred to another location to avoid over-writing of the data.
- Find on Potassium :
 - raw (alias)
 - `ls -lrt`
 - `cd Nov_15`
 - sub-directory tps
 - tps_00_06

Finding Data the Next Day on Potassium


- *All* data transfers from the IDBM Database to Potassium at 2am the next day *

Do NOT shutdown Signa and computers at the end of the day!

- Telnet to potassium:


```
telnet
open 128.231.212.156
login:sdn-nfs
password:xxxx
at the $ prompt, type
epirt
ls -lrt
cd DATA_Jun_30
    GINX (anat files) MROR (epiRT)
```

Help! Where is My Data 1?



- Pfiles:
 - Potassium
 - sdc-adw
- GE RT
 - Potassium
 - sdc-nfs
- Image Files (GINX)
 - SunAdw (GINX)
- Potassium
 - sdc-nfs

Help! Where is My Data 2 ?



- Image Files are saved to the MOD
 - These can be restored to Signa if deleted
- Transfer to the Sun ADW via the xferx (GINX) protocol:
 - Tar and send to Potassium (nfs)
- RT Folder (MROR)
 - {RT/2002.07.25/Ex6061rt.tar}
- Anat Folder (GINX)
 - {Anat/2002.07.25/Ex6061anat.tar}


Screen Lock Up on the Sun ADW

- Hold down the Stop on the left side of keyboard and hit the **a** key at the same time.
- Then type: boot and hit enter
- Wait for the login prompt
 - Login: sdc
 - Password: xxxxxx

Screen Lock Up on Potassium


- Restart the Potassium x-server by pressing:
 - cntrl + alt + backspace keys @ the same time
- Login: sdc-adw
- Password: xxxxxx
- Login: sdc-nfs
- Password: xxxxxx

Trouble Shooting the GE RT Sequence



- Normal message reminding scan operator to check prescan values
- Could be a software glitch

Trouble Shooting the NIH EPI Sequence



- Researcher forgot to change CV from :
 prep = 1
 To
 prep = 0
- Could be a real software error, check the message board

Summary 1

- Do NOT shutdown Signa, the ADW or Potassium at the end of the day.
- It is important to logout of Potassium when finished scanning.

Summary 2

From Signa

- GINX Image Data
 - Tar and store for 2 weeks on Potassium
 - RealTime Data
 - Tar and store for 2 weeks on Potassium
 - P Files are stored on Potassium for 1 week only!
- *** It is each researcher's responsibility to transfer their own data ***

Acknowledgements

3 Tesla Core Facility

NHLBI

NMR F

Duke List Serve for MR Technologists

GE Medical Systems


Avotec, Inc.

Medrad, Inc.

Berlex Laboratories

ACR Glossary of MR Terminology

*Thank You for Your
Attention*



National Institutes of Health