AFNI & SUMA

Concepts, Principles, Demos



















Released under the GNU General Public License Version 2 (GPL) [or any later GPL version]

AFNI is a research tool.

Clinical uses are not supported or advised.











http://afni.nimh.nih.gov/afni

Some Goals of FMRI Analyses

- Task-based experiments
 - Per subject: estimate amplitude of BOLD response to each different type of stimulus
 - Find+model inter-regional correlations
 between fluctuations in BOLD responses
- Resting-state experiments
 - Measure spatial patterns in coherent fluctuations in spontaneous BOLD
- Group level
 - Combine and contrast per subject results

pre-processing

- Time shifting = pretend get 3D snapshot
- Despiking = remove large blips
- Image Registration (AKA alignment)
 - intra-EPI time series, and EPI-Structural
- Blurring in space = lower resolution :-(& less noise :-) & more group overlap :-)
- Masking = ignore non-brain voxels
- Scaling = normalizing data amplitude
 - Makes inter-subject comparisons more valid

- Time series regression
 - model of the BOLD response in the data =
 Hemodynamic Response Function stimulus timing
 - plus baseline (null hypothesis) model
 - plus physiological noise
 - plus allowing for serial correlation
- Talairach-ing = Spatial Normalization
 - Talairach, MNI-152, ...
 - affine and nonlinear spatial transformations

- Group Analyses = Putting it all together
 - ANOVA, LME, Meta-Analyses, ...
- **Blobs** = Spatial models of activation
 - Assigning statistical significance to blobs
- Connectivity = Inter-regional analyses
 - SEM, PPI, SVAR, DCM, Granger, ...
 - Resting state FMRI (Connectome!)
- Dimensional factorization
 - Components, such as PCA, ICA, ...

- Data Formats = NIfTI-1.x is your friend
- Software for FMRI analyses: *ope
 - *open-source
 - AFNI*, BrainVoyager, FSL*, SPM*, ...
 - Whichever you use, don't blindly assume the software works perfectly all the time
- Most important thing I will say today
 Understand and check the steps
 applied to your data!
- 2nd most important: Is no "best" way to analyze data, just "reasonable" ways

AFNI = Analysis of Functional Neurolmages


- Developed to provide an environment for FMRI data analyses
 - And a platform for development of new software tools
- AFNI refers to both the program of that name and the entire package of external programs and plugins (more than 200)
- The Prime Directive in the development of AFNI:
 - Allow users to stay close to their data and view/analyze it in many different ways
- **SSCC** = Scientific Computing and Statistical Core
 - Our mission is help NIH (and beyond) investigators carry out the analyses of their (F)MRI data
 - Development of data analysis methods and putting them into usable and (reasonably) reliable software
 - Consulting and question answering and hand-holding





Fundamental AFNI Concept

- Basic unit of data in AFNI is the dataset
 - A collection of 1 or more 3D arrays of numbers
 - Each entry in the array is in a particular spatial location in a 3D grid (a voxel = 3D pixel)
 - Image datasets: each array holds a collection of slices from the scanner
 - Each number is the signal intensity for that particular voxel
 - Derived datasets: each number is computed from other dataset(s)
 - e.g., each voxel value is a t-statistic reporting "activation" significance from an FMRI time series dataset, for that voxel
 - Each 3D array in a dataset is called a sub-brick
 - There is one number in each voxel in each sub-brick.







Parts of **AFNI**

- Interactive visualization and analysis AFNI and SUMA
 - For looking at data and results
 - AFNI is based on 3D volumes = data as gathered by MRI
 - SUMA is based on 2D surfaces = models of cortical surfaces
 - A few kinds of analysis can be done by pointing+clicking
- Batch mode programs and scripts
 - Are run by typing commands directly to computer, or by putting commands into a text file (script) and later executing them
 - Most AFNI complex analyses are done in batch programs
- Plugins and Plugouts
 - Separate programs that attach themselves to AFNI and/or SUMA to provide extra capabilities



AFNI & SUMA Interlude





AFNI Batch Programs

- Many many important capabilities in AFNI are only available in batch programs
 - A few examples (of more than 100)
- Über-scripts: afni_proc.py and align_epi_anat.py
 - FMRI time series pre-processing and analysis
 - Driver for 3D image registration tools
- <u>3dDeconvolve</u> + <u>3dREMLfit</u> = multiple *linear* regression on 3D+time datasets; fits each voxel's time series to activation model, tests these fits for significance (<u>3dNLfim</u> = nonlinear fitting)
- <u>3dvolreg</u> = 3D+time dataset registration, to correct for small subject head movements, and for inter-day head positioning
- 3dANOVA + 3dLME + 3dMEMA = ANOVA/t-test group analyses: combining & contrasting datasets in Talairach space
- 3dsvm = SVM multi-voxel pattern analysis program
- 3dDWItoDT = compute diffusion tensor from DWI (nonlinearly)





Analysis by Super-Script – by hand

Script to analyze one imaging run (5 min) of data from one Subject [cd AFNI_data6/afni ; tcsh quick.s1.afni_proc]
 afni_proc.py -dsets epi_r1+orig -copy_anat anat+orig \ -tcat_remove_first_trs 2 \ -do_block align \ \

```
-execute
```


-regress stim times quick.rl times.txt

- Stimulus timing in file quick.r1_times.txt
 0 30 60 90 120 150 180 210 240 270
- 20 s of stimulus per block, starting at the given times
- FMRI data in file epi r1+orig
 - Anatomical volume in file anat+orig
- **Actions**: Align slices in time; align Anat to EPI; motion correct EPI; blur in space; activation analysis (thru time) in each voxel





Analysis by Super-Script – by GUI







FMRI Experiment Design and Analysis

FMRI experiment design

All on one unreadable slide!

- Event-related, block, hybrid event-block?
- How many types of stimuli? How many of each type? Timing (intra- & inter-stim)?
- Will experiment show what you are looking for? (Hint: bench tests)
- How many subjects do you need? (Hint: the answer does not have 1 digit)
- <u>Time series data analysis</u> (individual subjects)



- Assembly of images into AFNI datasets; Visual & automated checks for bad data
- Registration of time series images (AKA motion correction and EPI-anat alignment)
- Smoothing & masking of images; Baseline normalization; Censoring bad data
- Catenation of multiple imaging runs into one big dataset
- Fit statistical model of stimulus timing+hemodynamic response to time series data
 - o Fixed-shape or variable-shape response models [pattern matching in time]
- Segregation into differentially activated blobs (i.e., what got turned on or off?)
 - Threshold on statistic + clustering <u>and/or</u> Anatomically-defined ROI analysis
- Visual examination of maps and fitted time series for validity and meaning
- Group analysis (inter-subject)
 - Spatial normalization to Talairach-Tournoux atlas (or something like it; e.g., MNI)
 - Smoothing of fitted parameters
 - Automatic global smoothing + voxel-wise analysis or ROI averaging
 - ANOVA+ to combine and contrast activation magnitudes from the various subjects
 - Visual examination of results (usually followed by confusion)
 - Write paper, argue w/ mentor, submit paper, argue w/ referees, publish paper, ...





Getting and Installing AFNI

- AFNI runs on Unix systems: Linux, Sun, Mac OS X
 - Can run under Windows with Cygwin Unix emulator
 - This option is really just for trying it out not for regular use!
- If you are at the NIH: SSCC can install AFNI and update it on your system(s)
 - You must give us an account with ssh access
- You can download precompiled binaries from our Website
 - http://afni.nimh.nih.gov/afni
 - Also: documentation, message board, humor, data, ...
- You can download source code and compile it
- AFNI is updated fairly frequently, so it is important to update occasionally
 - We can't help you with old versions!





AFNI at the NIH Scanners

- **AFNI** can take 2D (or 3D) images in "realtime" from an external program and assemble them into 3D+time datasets slice-by-slice at each TR then update the images+graphs
- Jerzy Bodurka (ex-FMRIF) and Vinai Roopchansingh have set up the GE FMRI scanners (3 Ts, 1.5 T, and 7 T) to start AFNI automagically when scanning, and send reconstructed images over to the AFNI box as soon as they are available:
 - For immediate display (images and graphs of time series)
 - Plus: graphs of estimated subject head movement
 - Also possible: feedback to subject in realtime
- Goal is to let you see image data as they are acquired, so that if there are any big problems, you can fix them right away
 - Sample problem: someone typed in the imaging field-of-view (FOV) size wrong (240 cm instead of 24 cm), and so got garbage data, but only realized this too late (after scanning 8 subjects this way) D'oh!





That's All, Folks



